



Model CB1004

CARER BED

An electrically controlled 'pull-apart' Bed

FEATURES: This commercial quality 4-function bed is specifically designed for home care, rental and care facilities.

Ideal for both short term and extended care.

Comes in 4 easy to assemble parts, without the use of tools.

Robust in construction but light weight, with no part weighing more than 21kg

Head & foot boards allow good visibility and double as armrests to assist patients in & out of bed.

Easy storage - using the 'transporter frames' (optional), the whole bed can be stored vertically.

Compact height (1220mm) fits in most vans.

SPECIFICATIONS:

Mattress size:	900 x 1960mm
SWL	220 Kg
Height Range	345 - 745mm
Trendelenburg	0° - 11.3°
Wheels	100mm dual caster with total lock
Electric System	Dewert
Back rest	0° - 70°
Knee Break	0° - 45°
Paint Finish	Power Coated (White Birch)
Head & Foot board	Laminated Antique White
Total Weight	77.5Kg
Size in Transport Frame	570 x 931 x 1220mm



The Carer Bed has been designed & tested to meet AS/NZS: 3200.2.38:2007

ACCESSORIES:

Transporter Frame	Over bed table
Bed Stick	Self Help Pole & Socket
Side Rails	Remote control from phone app
Battery Back up	



Designed and manufactured by: **DEUTSCHER HEALTHCARE**

731 Creswick Rd • Ballarat 3350 • www.deutscherhealthcare.com.au • sales@deutscherhealthcare.com.au