

INTRODUCING A FOLDABLE FLOOR LEVEL BARIATRIC BED

with 350kg capacity



WMR1090BT (BIG TED) FOLDING BARIATRIC FLOORLINE BED

Welcome to the new generation in bariatric care, where you have all the advantages of a floor level bed with a 350kg SWL that folds. Installation & storage of this bed is so simple, added to this we offer side rails that can stay on the bed when folded.

Combined with all the features that Deutscher Health Care include as standard makes this bed the ideal package;

- Bracket for air mattress pump
- Built in bed control panel plus client hand set controller
- Extensive height adjustment from 135 to 740 mm
- Adjustable back rest, knee bend with twin lift actuators
- Removable head and foot boards coloured to match your décor or assist with asset management
- Large 125mm (5") wheels with total lock "easy use" control bars
- King single size mattress base of 1070 x 2030mm
- Mattress retainers

Australian made with 10-year warranty on the frame and 5 years on the electrics.

Affordable Quality built to both Australian and International standards. Full range of accessories, including - bed sticks, battery back-up, self-help pole, side rails and under bed lighting.



**NEED TO KNOW MORE,
TALK TO US WE ARE HERE
TO HELP**



DEUTSCHER HEALTHCARE

731 Creswick Rd • Ballarat 3350 • phone (03) 53395708 • sales@deutscherhealthcare.com.au