



comfort

DEUTSCHER HEALTHCARE

731 Creswick Rd • Ballarat 3350 • phone (03) 53395708 • sales@deutscherhealthcare.com.au

AUSTRALIAN MADE



BED STICKS

Bed sticks are utilised to assist the user to transfer to and from the bed and are also fitted to help a user move when in bed.

They can be easily fitted to either side of the bed without the use of tools.

These Bed sticks have a unique clamping design to ensure the bed stick is solidly mounted to the mattress base / bed frame. We manufacture a range of fixed pole or fold down pole bed sticks to suit most medical - aged care style beds.

Deutscher Clamp on bed sticks are designed for the user's safety and protection.

A 50mm ball on the end of the pole ensures the protection of the user in case of an accidental fall.

Our bed sticks are tested to withstand a 50kg load force in any direction.

Deutscher bed sticks are proudly manufactured at our facility in Ballarat, Victoria.

FOLDING

The WL series bed stick pole measures 510mm high from the top of the mattress base to the top of the 50mm ball. The WLR series measures 560mm and the WL..... T series bed sticks measure 590mm high.



WL1040F



WL1040FT



WL1045F



WLR1045F

BED STICKS

WL1040F / WL1040

Clamp is designed to suit 32mm x 50mm frames (common on high/low beds). The pole is mounted onto the side of the clamp allowing clearance for bed linen between the mattress and bed pole.

WL1040FT / WL1040T

As above however the pole is mounted onto the top of the clamp minimising the clearance between the mattress and the bed pole.

WL1045F / WL1045

Clamp is designed to suit 25mm x 50mm frames (common on high/low beds). The pole is mounted onto the side of the clamp allowing clearance for bed linen between the mattress and bed pole.

WLR1045F / WLR1045

Clamp is designed to suit 25mm x 25mm mattress frames common on floor line ultra-low beds.

WL1041 (FIXED BED STICK ONLY)

Clamp is designed to suit 40mm x 40mm frames found on some of the older style medical Hi – Low beds. The pole is mounted onto the side of the clamp allowing clearance for bed linen between the mattress and bed pole.

Clamp designs to suit most medical bed frame sizes.

FOLDING	MODEL	BED SIDE RAIL SIZE
	WL1040F	32mm X 50mm
	WL1040FT	32mm x 50mm
	WL1045F	25mm x 50mm
	WLR1045F	25mm x 25mm

FIXED POLE	MODEL	BED SIDE RAIL SIZE
	WL1040	32mm X 50mm
	WL1040T	32mm x 50mm
	WL1045	25mm x 50mm
	WLR1045	25mm x 25mm
WL1041	40mm x 40mm	

FIXED

The WL series fixed pole bed stick measures 560mm high from the top of the mattress base to the top of the 50mm ball. The WLR series also measures 560mm and the WL... ..T series fixed pole bed sticks measure 590mm high.



WL1040



WL1040T



This image illustrates both the WL1040F and the WL1040FT (foreground) mounted onto the bed frame. Note the distance between the bed stick pole and the mattress. The top mounted bed stick (T) reduces the mattress gap and is recommended for use if this gap is an issue



WLR1045



WL1045





731 Creswick Road
Ballarat, Victoria, 3350

T: 03 5339 5708

www.deutscherhealthcare.com.au

H F Deutscher Pty. Ltd. Have a policy of “maintaining Australian made content” and are proud to be an Australian owned and operated company producing products for over 110 years.