



comfort

DEUTSCHER HEALTHCARE

731 Creswick Rd • Ballarat 3350 • phone (03) 53395708 • sales@deutscherhealthcare.com.au

AUSTRALIAN MADE



MENTAL HEALTH BED - BATTERY WZ1005

The mental health bed has been specifically designed to meet the requirements of mental health institutions. It has a strong steel frame and welds that have been tested to safely support a rated load of 250kg. The lift actuator and electrical components are inset located and protected under the mattress base giving the bed a clean appearance. All wiring is completely contained in the steel work of the bed. The bed surface is a flat sheet with no backrest or knee bend options. The slim line head & footboards are not removable once installed. All critical moving joints are fitted with bronze bearings for quiet operation and long life. The height range of the bed is 240 – 740mm.

The bed is operated by the control mounted on the portable battery pack trolley, this unit is completely removable from the area once the bed has been adjusted.

The bed has a smooth tough powder coated finish, giving a durable surface and safe, clean appearance to the bed. This bed has been tested for over 2000 cycles with its rated load. Other testing has been carried out in accordance with Australian Standards to ensure that the bed is safe, reliable and meets the requirements of the TGA. Once again, being an Australian company, we are always only a phone call (or Email) away for any assistance that may be required.

- Options:
- Larger wheels
 - Head and Foot board colours

WZ1005	Single Bed	920x1960mm mattress base
WZ1006	Battery Pack	Includes charger and wall bracket
WZ1007	Trolley	Transport Trolley
WZ1008	Electric Pack	Trolley battery bracket, handset & connection cable