



comfort

# DEUTSCHER HEALTHCARE

731 Creswick Rd • Ballarat 3350 • phone (03) 53395708 • sales@deutscherhealthcare.com.au

AUSTRALIAN MADE



WM Series



DW Series



WM Series



## DEUTSCHER WALMSLEY RANGE

### WM Series

The Walmsley WM series is a commercial quality Hi Lo electric bed that is reliable and safe. With a maximum load capacity of 250kg, it is ideal for hospital, hostel and home care environments. It has an electrically operated back rest, knee bend, raise and lower functions, while the hand control has a self coiling lead and is easily cleaned.

The electrical system is located under the mattress base for safety and protection, while all critical moving parts are fitted with bronze bearings for quiet and smooth operation. The head & footboard are easily removed (without the use of tools) for safety, storage and transportation, while a tough powder coated finish provides a durable and fresh appearance that is easily cleaned and maintained.

Three mattress base sizes are available to suit client needs:

Standard size	920 x 1960mm
Standard wide	1040 x 1960mm
King size	1040 x 2030mm

The Walmsley beds have been designed tested and proudly manufactured in Ballarat, Australia to meet both Australian and International standards.

See chart on back page for specifications.

### DW Series

This bed has been designed mainly for home care users. It is reliable and trouble free but lacks some of the features & options of the WM series. The DW Series has an electrically operated backrest, knee bend & raise & lower, a generous mattress base (920 x 1960mm), and a maximum load capacity of 185kg.

The hand control is commercial quality, comes with a self coiling lead which is easily cleaned, while all critical moving parts are fitted with bronze bearings for quiet and smooth operation. The head board is easily removed for safety & storage, while a mattress retainer is provided at the foot end of the bed. (Foot board is available). A tough powder coated finish provides a durable and fresh appearance that is easily cleaned and maintained.

The DW series provides a great alternative for those who do not require the features of the WM series bed. This bed is proudly manufactured in Ballarat, Australia to meet both Australian and International standards.

See chart on the back page for specifications.

## Specifications

	DW Series (light weight)	WM Series (commecial)
Safe working load	185 KG	250KG
Wheels (std)	100mm Dual Brake & Lock	100mm Dual Brake & Lock
Wheels 125mm option	Yes	Yes
Central locking option	Yes	Yes
Electrical system	DEWERT	DEWERT
Battery backup option	N/A	Avaiilable
Lockout function on hand set	N/A	Avaiilable
Handset on flat easy to clean pad	STD	STD
Infra red handswitch	Avaiilable	Avaiilable
Underbed light	Avaiilable	Avaiilable
Leg raise	STD	STD
Raised backrest	STD	STD
Trendelenburg (TILT)	NO	NO
Height range of bed	24 -74cm	24 - 74cm
Hoist clearance	140mm	140mm
Removable head board	Yes	Yes
Removable foot board	OPTIONAL	Yes
Optional head & footboard designs	Yes	Yes
Rear corner buffers	Avaiilable	Avaiilable
Mattress flatform (std)	92 x 196cm	92 x 196cm`
Mattress flatform (std wide)	N/A	107 x196cm
Mattress flatform (king single)	N/A	107 x 203cm
Self help pole sockets	N/A	STD
Side rails - 3 or 4 rail	YES	YES
Alternating air mattress bracket	NO	NO
Paint finish powder coated std	White birch	White birch
Mattress retainers	N/A	STD
Warranty on frame	10yrs	10yrs.
Warranty on electrical	3yrs	3yrs
Weight of bed	70kg	80kg
Transporter frame	OPTIONAL	OPTIONAL
Fold for transportation	NO	NO
<b>Accessories:</b>		
Side Rails 3 bar	Avaiilable	Avaiilable
Side Rails 4 bar	Avaiilable	Avaiilable
Clamp on bed stick	Avaiilable	Avaiilable
Self help pole and handle	N/A	Avaiilable
Over Bed table (spring assisted height adjustment)	Avaiilable	Avaiilable
Foot extension	Avaiilable	Avaiilable
Bed side locker	Avaiilable	Avaiilable
Foot end mattress retainer instead of foot board	Std	Avaiilable



**WM1004**  
STANDARD SINGLE

**WM1014**  
WIDE SINGLE

**WM2040**  
KING SINGLE



**DW1004**  
STANDARD SINGLE

## ELECTRICALLY OPERATED SYSTEMS

Deutscher Healthcare Beds are electrically operated by Dewert Actuators. Dewert is a German electrical engineering company that is reputed to be the largest manufacturer of linear actuators in the world.

Dewert are well supported in Australia, with a large distribution warehouse in Melbourne and service agents operating throughout the country. All actuators used by Deutscher Healthcare are supported with test certificates, proving conformance with the relevant standards. Deutscher Healthcare is very selective in regard to the quality of systems used on its products. We research each product to ensure it is safe, fit for the purpose and able to carry out the tasks required reliably. All handsets used on our equipment are Keypad - these are easily cleaned and do not absorb moisture. All Dewert actuators used by Deutscher Healthcare are rated IPX4 or above for water resistance.

HF Deutscher Pty. Ltd. is an Australian, family owned and operated business, which has been manufacturing in this country for over 100 years. Since our inception in 1904 we've had a long and proud history of producing goods that have a positive effect on peoples day to day lives. In the early years we manufactured stationary engines and equipment for the farming sector, before diversifying into a range of commercial quality mowers - some of which are still manufactured today. In 1938 Henry Deutscher registered a patent for a mechanism used on a medical bed. Today, Deutscher Healthcare manufactures a range of homecare, aged care & medical beds. Deutscher Healthcare beds have been designed to comply with Australian Standards, ensuring that the bed you purchase is safe, reliable and suitable for your needs.

Deutscher Healthcare has a modern manufacturing facility, centrally located in Ballarat, Victoria with close proximity to major transport routes.

