



comfort

# DEUTSCHER HEALTHCARE

731 Creswick Rd • Ballarat 3350 • phone (03) 53395708 • sales@deutscherhealthcare.com.au

AUSTRALIAN MADE



MEETS THE DHHS REQUIREMENTS



## DEUTSCHER FLOORLINE LUCY BED

### WMR FLOORLINE SERIES

The new generation WMR Floorline bed is designed so the client can be placed low to the floor and raised to a height to allow carers to safely care for the client. The new slimline towers give the bed a clean, smooth and attractive appearance. This bed comes standard with adjustable backrest, knee bend, high & low height adjustment along with a tilt function, and can be easily folded for transportation, making it easy to be installed or stored. A control panel is located in the foot end tower of the bed, while a handset controller is offered as standard for patient control of the bed.

The WMR Floorline bed comes with a standard mattress size of 900 x 1930mm or in a king size version of 1070 x 2030mm for taller clients. A mounting support bar is fitted to the foot end of the bed, should an alternating mattress be required. The bed is rated to give a 220kg safe working load, while the transformer is built into the bed tower to prevent damage when lowered.

All our beds have been designed, tested and proudly manufactured in Ballarat, Australia to meet both Australian and International standards.

#### Order descriptions:

- WMR1080LC - Standard single size bed
- WMR1090LC - King single size bed

#### Options

- Bed Sticks
- Battery back-up
- Self help pole and adaptor socket
- Under bed light
- 50mm bolt in bed extender

- Maximum load rated to 220Kg
- Removable FULL width head & foot boards
- Mattress pump bracket
- Quality Dewert electrics
- Dual castor 100mm wheels
- Wheels fitted with brake & direction lock
- Mattress platform height range: from 125mm - 750mm + 150mm Mattress = 900mm (meets with DHHS requirements)
- Warranty: Frame 10 Years, Electrics 5 Years



*pictured above:  
Built-in Control Panel  
Left:  
side view of folded unit*

