



DEUTSCHER HEALTHCARE



FLOORLINE BED - ROSE II



WMR ROSE II

The WMR series bed allows patients to be lowered to the floor and raised to a suitable height when carers need to attend the client. It has electrically operated backrest, knee break and raise and lower and trendelenberg tilt.

Slimline towers at either end give the bed a clean, smooth and attractive appearance. The footend tower has a built in keypad to operate the bed. A separate hand control comes standard with the bed. This bed can be folded for storage, transport and installation and easily fits through any standard doorway.

The Rose II is available as a standard single 900 x 1930mm and king single 1070 x 2030mm mattress platform. It is designed and manufactured using quality materials in Ballarat Australia.





FLOORLINE - ROSE II

STANDARD FEATURES

WMR 1080 WMR 1090	Standard Single 900 x 1930mm King Single 1070 x 2030mm	Mattress Base sizes
Height Range	125 to 750mm (Top of mattress platform)	
Wheels	Dual castor brake and direction lock 100mm	
Pump Bracket	Fitted on the footend tower	
Finish	Powder coating – Pearl White	
Safe Work Load	220kg (includes client and bedding)	
Controls	Handset 5499 and keypad built into footend tower	
Warranty	10 years for Frame 3 years for electrical system	
Head and Foot Boards	Standard Cherry Red. Other colours available	
Unfolded length	Single 2300mm King Single 2400mm	

OPTIONS AVAILABLE

WLR1050 S	3 bar side rails Single Bed
WLR1055 S	4 Bar side rail Single Bed
WLR1050 K	3 bar side rails King Single Bed
WLR1055 K	4 bar side rails King Single Bed
WLR1060	Self help pole and handle
WLR1062	Socket
WLR1045	Straight bed stick
5869	Mattress retainers
WLR1045 F	Folding bed sticks
Light	Under bed lighting
WMR BAT	Battery backup
5827	CPR Function and magnetic lock out
5866	Back light and magnetic lock out
WMR1071	50mm Bolt in Mattress Base Extender

July 2017