


Deutscher Healthcare

Over bed tables and bed side lockers-

Over bed table #WL1030



- Table top size 750 x 400mm
- Lower trolley inside size 690 x 390mm
- From floor to top of lower trolley 100mm
- Floor to table top-
 - Minimum 740mm
 - Maximum 1050mm
- SWL 20 Kg in centre of table
- Over all weight 11 Kg
- Wheel diameter 50mm

Bed side locker #5686 3 draw



3 side return safety table top

- Table top size W425 x L450mm Inside
- Over all dimensions W465 x L475 x H765mm (H735 to top)
- Inside top draw dimensions W375 x L375 x D140mm
- Internal second draw dimensions W375 x L375 x D140mm
- Internal lower draw dimensions W375 x L375 x D215mm
- Wheel diameter 60mm
- SWL on table top 25Kg
- SWL top draws 5Kg (Opening one draw at a time)
- SWL lower draw 15Kg (extra safety wheel under lower draw)
- Over all weight 34Kg

Bed side locker #5620 4 draw



Flat table top

- Table top size W505 x L480mm
- Over all dimensions W460 x L450 x H840mm
- Inside draw dimensions W375 x L375 x D140mm
- Wheel diameter 60mm
- SWL on table top 40Kg
- SWL draws 5Kg (opening one draw at a time)
- Over all weight 39Kg