



**DEUTSCHER  
HEALTHCARE**  
*Manufacturing Engineers since 1904*

## High Protection System (HPS) for Deutscher Floor Line Bed



**High Protection System for the Deutscher floor-line bed that offers maximum comfort and protection- WLR1080HPS.**

The High Protection System (HPS) is designed to offer maximum comfort and protection to any patient suffering from continuous movement, including Huntington's, Cerebral Palsy, Acquired Brain Injury, Psychosis and Epilepsy.



- Strong but soft padded side panels and bed ends for user protection
- Hinged fold down side panels allow greater access to the user
- The sides can be completely removed, if required for even greater access when utilising patient hoists or when better access to the user is required
- The complete High Protection System can be removed when not required, converting the bed back to a standard floor level bed

## CLINICAL AND TECHNICAL DETAILS

- **Standard single size bed**
- **Lowest mattress base platform position** : 125mm
- **Highest mattress base platform position** : 750mm
- **Mattress Length** : 1930mm
- **Mattress width** : 900mm
- **Lifting / SWL Capacity:** 220Kgs  
(ref max patient weight 160kg)
- **Side panel height from mattress base** : 440mm
- **Head and foot board height from mattress base** : 590mm
- **Wheels** : 100mm diameter, twin wheel, individual tab control total locking casters
- **Overall Bed Weight:** 128Kgs
- **External Bed Dimensions:** 2250mm x 1190mm  
(1300mm wide including safety hooks)

For Further Information please contact:



Richard Jeffery

M: 04386 04106

P: (03) 9191 6517

E: [richardj@smtrbiocare.com](mailto:richardj@smtrbiocare.com)

W: [www.smtrbiocare.com](http://www.smtrbiocare.com)



Deutscher Healthcare

M: 0417 369964

P: (03) 53395708

E: [sales@deutscherhealthcare.com.au](mailto:sales@deutscherhealthcare.com.au)

W: [www.deutscherhealthcare.com.au](http://www.deutscherhealthcare.com.au)