



DEUTSCHER HEALTHCARE



WZ1005

MENTAL HEALTH BEDS

The Australian designed and manufactured Deutscher Healthcare MENTAL HEALTHCARE bed range has been developed to meet the requirements for use in all aspects of Mental Health Care including use in acute and intensive extended care psychiatric facilities.

The beds are available in three models:

- WZ1005** Single function Hi-Low height adjustable bed
- WZ1000** Two function Hi-Low height adjustable bed with adjustable back rest
- WZ1004** Three function Hi-Low height adjustable bed with adjustable back rest and knee bend

The **WZ1005** is the most popular bed for use in psychiatric care facilities as it is only single function (height adjustment) and low ligature risk with a completely removable power lead / power supply system utilising a battery pack or mains electricity adaptor to power the bed.

All beds have fixed head and foot boards that are only removable with the use of a tool, minimal LOW VOLTAGE wiring that is predominately encased in metal tube conduit and the beds are of a robust and safe design.

All beds can be optioned to suit your exact requirements including selections like: central locking wheels, built in attendant control / operation panel, handset control, battery power or mains power operation and even colour choices of head and foot boards.

We welcome the opportunity to speak with you regarding your bed requirements.





MENTAL HEALTH BED RANGE



WZ1005



WZ1000

With optional built in attendant control panel



WZ1006
Battery and charger station
• Battery • Wall mount Battery holder
• Charger



WZ1007
Battery pack trolley



WZ1008
Battery holder (no battery)
for trolley mounting
• Battery holder cradle • Bed
controller • Bed connection lead



WZ1010
Complete battery kit for trolley mounting
• Battery • Battery cradle (holder)
• Charger • Bed Controller
• Bed connection lead



WZ1009
Mains power adaptor
and control unit
(In place of battery)

SPECIFICATIONS

Mattress Platform: 196 x 92cm (STD Single bed)

Safe working load: 220kg. This includes patient weight & accessories

Finish: Powder Coated in Pearl White

Wheels: All wheels are 100mm diameter dual castor with independent tab operated wheel brake locks

Height Range: 240 - 740mm.

Electrics: DEWERT (Germany)

ABOUT DEUTSCHER HEALTHCARE

Deutscher Healthcare beds & accessories are manufactured in Ballarat, Victoria. Deutscher only use quality materials to ensure you are purchasing a quality bed that is designed to give trouble free service.

Deutscher have a modern facility located close to major transport routes.

Our sales staff are always keen to assist you with any further information you may require. The company has been in operation since 1904 and is still owned and operated by the same family.

All products are designed and manufactured to a high quality using high grade materials to ensure the product you buy meets your requirements. We ensure the bed you purchase is safe reliable and designed to give trouble free service for many years.

September 2019

731 Creswick Road, Ballarat, 3350 | (03) 53395708 | deutscherhealthcare.com.au

