



DEUTSCHER HEALTHCARE

HIGH PROTECTION SYSTEM (HPS) FOR DEUTSCHER FLOORLINE BED

WMR1080HPS-Single / WMR1090HPS-King Single

NEW WMR1090BTHPS-King Single "BT" Bariatric



***NEW**
NOW AVAILABLE
IN A BARIATRIC
VERSION
350Kg
SWL bed.



High Protection System (HPS)

for the Deutscher floorline bed offers maximum safety, protection & comfort.

The **High Protection System** (HPS) is designed to offer maximum protection to any client suffering with the effects of Huntington's, Cerebral Palsy, Acquired Brain Injury, Psychosis and Epilepsy.

Features:

- Strong but soft padded side panels and bed ends for user protection
- Hinged fold down side panels allow greater access to the user
- The sides can be completely removed, if required for even greater access when utilising client hoists or when better access to the user is required
- The complete **High Protection System** can be removed when not required, converting the bed back to a standard floor level bed .

CLINICAL AND TECHNICAL DETAILS

Lowest mattress base platform position	125mm / BT - 145mm
Highest mattress base platform position	750mm / BT - 770mm
Mattress Length	S - 1930mm / KS & BT - 2030mm
Mattress width	S - 900mm / KS & BT - 1070mm
Lifting / SWL Capacity	220Kgs (ref max patient weight 160kg) BT - 350Kgs (Max patient weight 260Kg)
Side panel height from mattress base	440mm
Head & foot board height from mattress base	590mm
Wheels	100mm diameter, twin wheel, individual tab control total locking castors BT - 125mm diameter, twin wheels with semi central locking system
Overall Bed Weight	S - 128Kgs / KS - 140Kgs / BT 160Kgs
External Bed Dimensions	S - 2250mm x 1300mm KS - 2350mm x 1470mm
Dual Controls	Handset and carer control panel



ABOUT DEUTSCHER HEALTHCARE

Deutscher Healthcare is a manufacturing company based in Ballarat, Victoria and has been in continuous operation for over 114 years. In earlier years, the parent company, HF. Deutscher Pty. Ltd., specialised in farm equipment and later moved into producing a range of mowers and outdoor power equipment, some of which is still manufactured today. The company is still owned and operated by the original family.

All products are designed and manufactured to a high quality, using first grade materials to ensure the product you buy meets Australian standard requirements; is safe, reliable and built to give trouble free service for many years to come.

For more information, speak with one of Deutscher's friendly staff about your needs and requirements.

Or, scan the QR code to visit our website now.



deutscherhealthcare.com.au



(03) 5339 5708



JULY 2023