



# DEUTSCHER HEALTHCARE



## BED STICKS

- Bed Sticks are used to assist the user to transfer in and out of the bed. They can also give assistance with in-bed mobility.
- Deutscher Bed sticks are designed to accommodate a wide range of medical beds used throughout the Healthcare industry.
- This range includes fixed and folding bed sticks.
- A 50mm ball fitted to the top of the stick protects the user should they fall and provides a good gripping surface for the client.
- All Deutscher bed sticks are securely clamped to the mattress platform.
- They can be fitted to either side of the bed.
- All Deutscher bed sticks are tested to withstand a 50kg force in all directions ensuring it is safe to use.
- All Deutscher bed sticks are designed and manufactured in Ballarat, Australia.



# BEDSTICKS



WL1040 & WL1040T  
WL1045 & WL1045T



WL1044 & WL1041




WL1044H, WL1045H  
& WL1040H



WLR1045 & WLR1045F  
WLR1046 & WLR1046F



WL1040FT & WL1045FT



WL1040F & WL1045F

PART/NO	FIXED POLE BED STICK RANGE	RAIL SIZE ON BED	SUITABLE FOR
WL1040	Side Mounted Fixed Pole bed stick	32 x 51mm	Walmsley Bed
WL1040T	Top Mounted Fixed Pole Bed Stick	32 x 51mm	Walmsley Bed
WL1045	Side Mounted Fixed Pole Bed Stick	25 x 50mm	Carer Bed
WL1041	Side Mounted Fixed Pole Bed Stick	40 x 40mm	
WLR1045	Side Mounted Fixed Pole Bed Stick	25 x 26mm	Floorline Bed
WL1044	Side Mounted Fixed Pole Bed Stick	20 x 50mm	Etude Bed
WL1044H	Side Mounted Fixed Pole Bed Stick (Hook Top)	20 x 50mm	Etude Bed
WL1040H	Side Mounted Fixed Pole Bed Stick (Hook Top)	32 x 51mm	Walmsley Bed
WL1045H	Side Mounted Fixed Pole Bed Stick (Hook Top)	25 x 50mm	Carer Bed
WL1040LX	Special Top Mounted Fixed Pole Bed Stick	32 x 51mm	Lexia Bed
WL1040LXH	Special Top Mounted Fixed Pole Bed Stick (Hook Top)	32 x 51mm	Lexia Bed
WLR1046	Side Mounted Fixed Pole Bed Stick	28 x 28mm	Floorline (Plastic Mattress sheet)
WLR1045H	Side Mounted Fixed Pole Bed Stick (Hook Top)	25 x 26mm	Floorline

PART/NO	FOLDING POLE BED STICK RANGE	RAIL SIZE ON BED	SUITABLE FOR
WL1040F	Side Mounted Folding Pole Bed Stick	32 x 51mm	Walmsley Bed
WL1040FT	Top Mounted Folding Pole Bed Stick	32 x 51mm	Walmsley Bed
WLR1045F	Side Mounted Folding Pole Bed Stick	25 x 26mm	Floorline Bed
WLR1046F	Side Mounted Folding Pole Bed Stick	28 x 28mm	Floorline (Plastic mattress Sheet)
WL1045F	Side Mounted Folding Pole Bed Stick	25 X 50mm	Carer Bed
WL1040LXF	Special Top Mounted Folding Pole Bed Stick	32 x 51mm	Lexia Bed

Top mounted bed sticks are often used in conjunction with side rails or where the bed stick needs to be firmer against the matterress.